



# vestige

## IMPLA CLEAR

Vestige IMPLA CLEAR is a transparent medium consistency addition silicone, suitable for direct impressions in the oral cavity and for the creation of transparent masks.

### Description / Use

Vinylpolysiloxane for the making of transparent masks, occlusal key, bite rebasing, and model stabilisation in guided implantology.

#### Standard packaging REF 1S6100

VESTIGE IMPLA CLEAR: 2x50ml cartridges + 6 mixing tips

### Characteristics

Fluid  
Transparent  
Rigid after setting  
Dimensional stability  
Safe to disinfect  
Easy to mix

### Advantages

Flowing, it enables reproduction without incorporating bubbles  
It does not interfere with resin polymerization and reduces its inhibited layer  
It does not deform  
Accurate and unaltered in time impressions  
Safe for dentists  
Packaged in standard 50ml cartridge

TECHNICAL DATA	VESTIGE IMPLA CLEAR
ISO 4823 / ADA spec. 19	Type 2 medium bodied consistency
Mixing ratio	1:1
Mixing time (23°C)	Automix
Total working time (23°C)	40"
Time in mouth (35°C)	2'
Setting time (35°C)	3'
Hardness	70 ± 2 Shore "A"
Linear dimensional change	< 0,20%
Elastic recovery	> 99,5%