

VESTIGE LAB FLOW ADVANCED

**VESTIGE LAB FLOW ADVANCED IS AN EXTREMELY PRECISE ADDITION SILICONE.
IT IS A THIXOTROPIC CORRECTION FLUID CREATED TO BE USED IN PROSTHODONTIC LABORATORY.
NORMAL SETTING TIME, EXTREMELY PRECISE IMPRESSION RECOGNITION ON THE MODEL.**

DESCRIPTION

Medium consistency vinylpolysiloxane destined for applications that require moulding with mask rebasing. Normal working/setting time. It can be scanned. Color: purple.

PACKAGING

Standard packaging REF 1L3000

VESTIGE LAB FLOW ADVANCED: 2x50ml cartridges + 6 mixing tips

TECHNICAL DATA	VESTIGE LAB FLOW ADVANCED
Mixing ratio	1:1
Mixing time (23°C)	Automix
Total working time (23°C)	2"
Setting time(23°C)	6"
Hardness	56
Linear dimensional change	< 0,20%
Elastic recovery	> 99,5%

Limitation of the responsibility The indications related to the products TRAYART, both in oral form and written, are founded on the knowledge and experience acquired in dental field. The professional users are kept however to follow the indications for the use and to verify the fitness of the product according to the type of application. The professional user is assumed all the responsibilities for the consequential risks from a non correct use of the products TRAYART. A possible defectiveness of the products limits the responsibility of the manufacturer to the value of the products used by the professional.

FOR DENTAL USE ONLY